

## Big Case Back Table

# Renew!

program

The Big Case Back Tables are extremely well constructed and built for years of reliable service but general maintenance and reconditioning of table components should be considered to keep your table in optimal condition.

OR Specific's **Refurbish Program** gives you an opportunity to assess the condition your Big Case Back Table and request one of our sales representatives to update and restore your table.

### Refurbish Kit includes:

- **New** Four Casters (2 standard and 2 brake)
- **New** Clevis Pins for Shelf (2 each)
- **New** "T" handles for shelves (2 each)
- **New** 3 Heavy-Duty Drapes
- **New** Velcro (pre-cut replacement stripes)

429-RK Refurbish Kit for 6' table

428-RK Refurbish Kit for 5' table

427-RK Refurbish Kit for 4' table

### **NEW** Table Sizes Available

To meet a variety of operating room needs and configurations, OR Specific has engineered and now offers tables of 6ft, 5ft and 4ft lengths.



427 Big Case Back Table  
4ft Standard